Name:

Name:



- Test each substance in the table with red cabbage indicator:
 - o Put a small sample of the substance in a dimple
 - o Add a couple of drops of indicator
 - o Record the colour and pH in the results table

	Tooth- paste	Fizzy Juice	Bicarb of soda	Lemon Juice	Lime Juice	Liquid Soap	Yoghurt	Milk	Water
	3	BIOD	***************************************					milk	V W
Colour									
рН									

Red Cabbage Indicator



- Test each substance in the table with red cabbage indicator:
 - o Put a small sample of the substance in a dimple
 - o Add a couple of drops of indicator
 - o Record the colour and pH in the results table

	Tooth- paste	Fizzy Juice	Bicarb of soda	Lemon Juice	Lime Juice	Liquid Soap	Yoghurt	Milk	Water
	3	1 8100	* *					milk	V a
Colour									
рН									