
- Test each substance in the table with universal indicator:
 - o Put a small sample of the substance in a dimple
 - o Add a couple of drops of indicator
 - o Record the colour and pH in the results table



	Tooth- paste	Fizzy Juice	Bicarb of soda	Lemon Juice	Lime Juice	Liquid Soap	Yoghurt	Milk	Water
		BIOD	***					milk	9
Colour									
рН									

Universal Indicator

- Name:
- Test each substance in the table with universal indicator:
 - o Put a small sample of the substance in a dimple
 - o Add a couple of drops of indicator
 - o Record the colour and pH in the results table

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
acid-	←					n	eutr	al		_				→	alka

	Tooth- paste	Fizzy Juice	Bicarb of soda	Lemon Juice	Lime Juice	Liquid Soap	Yoghurt	Milk	Water
	3	BIOD	***************************************					milk	• •
Colour									
рН									